

Female Catheterization Simulator LF00856U INSTRUCTION MANUAL





ABOUT THE SIMULATOR

The **Life/form**® Female Catheterization Simulator is designed to duplicate the human condition as closely as plastics technology allows. Care and treatment should be used since abuse or rough treatment will damage the simulator almost as it would injure or cause pain to a patient. Note: To avoid the possibility of leakage, make sure to use a 16 French Foley catheter (supplied with the simulator).

LIST OF COMPONENTS

- A. Life/form® Female Catheterization Simulator
- B. Fluid Administration Set
- C. Plastic Positioning Pad
- D. 16 French Foley Silicone Catheter
- E. Life/form_® Simulator Lubricant, 2 oz.
- F. Small White Towelette
- **G.** Obstruction Tool
- H. Hard Carry Case




Figure 1

Figure 2

SET UP

A. To fill the Simulated Bladder

- Hang the fluid administration set 18" above the simulator. Use Fluid Supply Stand (LF01022U shown, sold separately) or similar device.
 Note: Hanging the fluid administration set higher increases the water pressure and may cause leaking.
- 2. Close the slide clamp on the tube of the fluid administration set.
- **3.** Squeeze open the top of the bag and fill with water. (To extend the useful life of your simulator, Nasco recommends distilled water.)
- **4.** Connect the quick disconnect fittings of the administration set and simulator tubing. *(See figure 1.)*
- **5.** Open the slide clamp on the administration unit and the pinch clamp on the simulator drain tube.
- 6. Allow water to flow through the simulator until it exits the drain tube. Close pinch clamp. (See figure 2.) The system is now filled with water and pressurized.

Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.

B. Catheter Selection

- Use only 16 French Foley catheters. The Female Catheterization Simulator is precisely
 engineered to provide the most lifelike experience possible. Use of the 16 French Foley
 catheter facilitates catheterization and helps minimize damage to the simulated urethra.
 Use of larger catheters results in unnecessarily difficult catheterization procedures. Use of
 smaller sizes will result in leakage around the catheter.
- 2. Only silicone catheters should be used. Latex catheters, if left in contact with the simulator, will damage the simulator and void the warranty.

C. Lubrication

- 1. Lubricate the simulator's urethra and the catheter liberally EVERY TIME a catheterization procedure is attempted. This provides maximum realism and minimizes damage to the urethra wall of the simulator.
- 2. Use only *Life/form*[®] Simulator Lubricant (LF00985U). Make NO OTHER substitutions. Even water-soluble lubricants tend to build up in the urethra over a period of time.
- **3.** Initial use after long storage or heavy usage during a long class session will require heavy lubrication. Do this by inserting the lubricated catheter about half way, then remove and relubricate several times before inserting the catheter fully.




Figure 3

Figure 4

GENERAL INSTRUCTIONS FOR USE

Catheterization Procedure

Note: To avoid the possibility of leakage, be sure to use a 16 French silicone Foley catheter (supplied with simulator).

- **1.** Position the plastic pad between the legs and under the buttocks to allow the proper angle for catheterization.
- Lubricate catheter by coating with Life/form® simulator lubricant. (See figure 3.) (Obtain a catheter bag or water-safe container; not included.)
- Slowly insert the catheter through simulated sphincter. With experience, you will feel a "pop" as the catheter passes into the bladder. The simulator will require about 3½" of catheter. Continue insertion until water flows from catheter. (See figure 4.)
- 4. Remove the catheter from the urethra after use.

Note: Special care should be taken when using a Foley catheter. Cuff inflation should only be attempted when the catheter is in the proper position inside the bladder. Just as in a real patient, the cuff must be completely deflated before the catheter is removed. Improper use of a Foley catheter can result in damage to the simulator.

TROUBLESHOOTING

- Lack of sufficient lubrication may make inserting the catheter difficult. See Set Up C instructions.
- Water does not flow even after proper insertion. Check fluid administration set water level. If there is still no water, remove catheter and repressurize the system to fill the Simulated Bladder. See Set Up A instructions
- **Dried Lubricant Obstruction.** Use the supplied obstruction tool to clear dried lubricant from the simulated urethra. Follow Lubrication and Catheterization Procedures using the obstruction tool. Flush the system with warm water to completely clear the obstruction and any added lubrication.

CARE AND MAINTENANCE

A. Cleaning the Simulator

- 1. Normal soil accumulated on the surface of the simulator can be removed with mild soap and lukewarm water.
- 2. Use Nasco Cleaner (LF09919U, sold separately) to remove stubborn stains from simulator. Simply apply to soiled area and wipe clean with a soft cloth or paper towels.

B. Storing the Simulator

Do not store simulator for extended periods of time with water in the unit. Before returning simulator to the case, drain as much liquid as possible from the bladder.

- 1. Open both tubing clamps. Drain the simulator by holding it upside down and high in the air.
- 2. Allow the water to drain out from the simulator.
- **3.** Dry all surfaces to prevent mildew formation while stored in the case.

CAUTION

Solvents or corrosive materials will damage the simulator. Never place simulator on any kind of printed paper or plastic. These materials will transfer indelible stains. Ballpoint pens will also make indelible stains.

SUPPLIES AND REPLACEMENT PARTS

LF00844U Obstruction Tool Kit (1 obstruction tool and 2 oz. lubricant)

LF00985U Life/form® Simulator Lubricant Kit

- LF01022U Fluid Supply Stand
- LF01098U Fluid Administration Set


LF01127U Foley Urethral Catheters (pkg. of 1)

- LF01128U Foley Urethral Catheters (pkg. of 10)
- LF09919U Nasco Cleaner

Other Available *Life/form* Simulators

LF00698U	Adult Injectable Arm (Light)	LF01038U
LF00855U	Male Catheterization	LF01062U
LF00898U	Ostomy Care Training	LF01068U
	Models Set	
LF00901U	Prostate Examination	LF01069U
LF00906U	Ostomy Care	LF01070U
LF00957U	Enema Administration	LF01084U
LF00980U	Advanced Breast Exam	LF01143U
LF01029U	Cath-Ed 1 and Cath Ed 2 Set	LF01193U
LF01035U	Infant Male & Female	
	Catheter Trainer	

LF01038U	Episiotomy Suturing Set	
LF01062U Pelvic, Normal & Abnormal		
LF01068U	Cervical Dilatation and	
	Effacement Simulators	
LF01069U	Cervical Effacement	
LF01070U	Birthing Station	
LF01084U	LF01084U Sigmoidoscopic Examination	
LF01143U	Testicular Exam	
LF01193U	Special Needs Baby	







ONLINE: EMAIL: PHONE: ADDRESS: NascoHealthcare.com lifeform@enasco.com 1.800.431.4310 901 Janesville Ave. Fort Atkinson, WI 53538