




Hemodialysis Practice Arm LF01037U Instruction Manual



WARNING: Product contains dry natural rubber. Do not remove film from tubing!

Life form Products by Nasco

About the Simulator

This Life/form® Hemodialysis Practice Arm Simulator duplicates a human arm. This special simulator contains an arteriovenous fistula. All aspects of the simulator are as lifelike in every detail as modern plastics technology allows. Treatment of the simulator must also duplicate that of a human arm. Rough handling will lead to damage.

This simulator is designed for student and patient teaching/demonstration. Its importance is in facilitating the teaching of techniques of easy access to the arm's blood supply for home dialysis.

It is possible to replace the skin and veins when it becomes necessary. The skin can be easily removed and new veins can be easily reinserted, providing you with a new arm. The life of the restored arm will be prolonged by utilizing only the size needles indicated in the instructions. The Skin and Vein Replacement Kits are available (see list of supplies).

Internal Structure

Internally, the vascular structure (rubber tubing) for the arteriovenous fistula begins in the inner aspect of the upper arm, travels down the arm, makes a loop in the inner aspect of the lower forearm, and then returns to the upper arm.

List of Components

- 2 IV bags
- Small towel
- Container of Life/form_® Hemodialysis Arterial Blood (LF01048U)
- · 10 cc syringe with needle
- · 3 cc syringe with needle
- Extra needle

Preparing the Artificial Arterial Blood and Vascular System

- Add 473 cc of distilled water to the container labeled Hemodialysis Arterial Blood.
- Pour the Hemodialysis Arterial Blood into one of the IV bags. Hang the bag no more than 18" above the level of the arm.
- 3. Attach the end of the IV bag to the arm tubing (the end without the pinch clamp).
- 4. With the other tubing in an appropriate receptacle, open the clamp on the IV tubing and the pinch clamp, and "flush" the system with the Hemodialysis Arterial Blood. Make sure all air bubbles have been removed from the system.
- 5. Once the "blood" has passed through the entire system, use the pinch clamp to close off the blood outlet tubing. Make sure the clamp on the IV tubing is open at this time. The system is now full of arterial "blood" and pressurized.
- The arm is now ready for demonstration and/or practice of easy access to the arterial system.
 The standard procedure for venipuncture can be utilized at this point. Insert only small diameter needles (20-25 gauge). This will prolong the life of the skin and veins. (In a patient, needles ranging from 15-17 gauge would be used).

Care of Simulator

After each class, disconnect the "blood" and flush the system. Unclamp the pinch clamp at the end of the tubing and remove IV sets from the arm. Use tap water to flush the system. Wash outside of arm with mild soap and water. Excess water may be removed from the fistula by raising the hand, lowering the upper arm, and draining it into a receptacle. Always drain the excess water from vascular system before storing.

Normal soil accumulated on the surface of the simulator can be removed with mild soap and lukewarm water. Use Nasco Cleaner (LF09919U) to remove stubborn stains from the simulator. Simply spray soiled area and wipe clean with a soft cloth or paper towels.

Solvents or corrosive materials will damage the simulator. Never place on any kind of printed paper or plastic. These materials will transfer indelible stains. Ball-point pens will also make indelible stains.




Skin Replacement

A. Remove Skin

- 1. Until lace from base of arm.
- 2. Lubricate outside of skin with Life/form_® Lubricant Kit (LF00985U).

- 3. Peel off skin carefully, turning inside out.
- 4. Remove all tubing.



B. Install New Tubing

1. Place tubing in channel of core.

C. Lubricate New Skin

 Pour lubricant into skin and swirl around so it covers all surfaces.





D. Slide Skin Over Hand of Core

- 1. Make certain the palm of the hand and core are in the same position.
- E. Grasp the Skin with Both Hands and Slide Over the Core Until the Fingers of Core Fit Into the Finger Holes





F. Work Fingers in Place



G. Check Tubing Position

 If tubing slipped from channel, it can usually be pushed back into place by working it with the fingertips from outside of skin.



Cautions

- The Life/form® Hemodialysis
 Arterial Blood is specially formulated to be compatible with the
 special veins and plastics used
 in manufacturing the arm.
- DO NOT use dull or burred needles, as these will cause leaks in the system.
- 3. DO NOT allow "blood" to dry on simulator or it may stain arm.
- Use only 500 cc of infusion fluid, as a larger amount will also increase the pressure of the venous system, resulting in leaks.
- DO NOT clean the simulator with solvents or corrosive materials, as they will damage it.
- 6. DO NOT use for subcutaneous or intramuscular injections. The Life/form® Intradermal Injection Simulator (LF01008U) and Adult Venipuncture and Injection Training Arm (LF00698U and LF00997U) are specially designed for this.

7. The Vein Tubing Sealant Kit (LF01099U) will extend the useful life of the tubing.

Supplies/Replacement Parts for Hemodialysis Practice Arm

LF01048U Life/form_® Hemodialysis

Arterial Blood — 1 Quart

Skin and Vein Replacement

LF01099U Vein Tubing Sealant Kit

Kit

LF01040U

LF09919U Nasco Cleaner

LF00985U Life/form® Lubricant Kit

Other Available *Life/form*. Simulators

Other Available <i>Life/form</i> 。Simulators			
LF00698U	Adult Injectable Arm (White)	LF01184U	Venatech IM & Sub Q
	Male Catheterization	LF01193U	Special Needs Baby
LF00856U	Female Catheterization		CPARLENE® Series
LF00901U	Prostate Examination	LF03601U	Adult Airway Management
LF00906U	Ostomy Care		Trainer with Stand
	Surgical Bandaging	LF03602U	Adult Airway Management
	Enema Administration		Manikin
LF00958U	Pediatric Injectable Arm	LF03609U	Child Airway Management
	Intramuscular Injection		Trainer with Stand
	Breast Examination		Child CRISIs™ Manikin
LF00995U	Arterial Puncture Arm	LF03617U	Deluxe Child CRISis ™
LF00999U	Pediatric Injectable Head		Manikin with Arrhythmia Tutor
LF01005U	First Aid Arm		PALS Update Kit
LF01008U	Intradermal Injection Arm	LF03623U	Infant Airway Management
LF01012U	Heart Catheterization (TPN)	1502/2011	Trainer with Stand
LF01019U	Ear Examination	LFU3032U	Child Intraosseous Infusion/ Femoral Access Leg on a Stand
LF01027U	Peritoneal Dialysis	I E0343311	Child Airway Management
LF01028U	Suture Practice Arm	LF03033 0	Trainer Torso
LF01034U	Suture Practice Leg	I F03693U	Basic Buddy® CPR Manikin
	Spinal Injection		"Airway Larry" Airway
	Hemodialysis Practice Arm		Management Trainer
	Episiotomy Suturing Set	LF03709U	Infant CRISis™ Manikin
LF01042U			Baby Buddy™ Infant CPR Manikin
	Pelvic, Normal & Abnormal		Fat Old Fred
	Stump Bandaging, Upper		Airway Management/Cricoid
	Stump Bandaging, Lower		Pressure Trainer
	Cervical Effacement	LF03770U	Chest Tube
	Birthing Station		CRiSis™ Manikin, Complete
	Cricothyrotomy	LF03955U	Deluxe <i>CRiSis</i> ™ Manikin
	Tracheostomy Care		Deluxe "Plus" CRiSis™ Manikin
LF01084U	Sigmoidoscopic	LF03965U	Adult CRiSis™ Auscultation
150100711	Examination		Manikin
	Central Venous Cannulation	LF03966U	Adult CRISIS™ Auscultation
	Blood Pressure Arm		Manikin with ECG Simulator
	Infant Intraosseous Infusion		GERi™/KERi™ Manikin Series
	Advanced IV Arm	LF04200U	Adult Sternal Intraosseous
	Venipuncture and Injection Arm	150/00111	Infusion
	Advanced IV Hand Auscultation Trainer	LFU0UU I U	CPR Prompt® Adult/Child Manikin
	Testicular Exam	I E0601211	CPR Prompt® Infant Manikin
	Male & Female Catheter	I FOA 20012U	CPR Prompt® Keychain
	Advanced CPR Dog	LFUUZUUU	Rescue Aid
	Venatech IV Trainer	LE0620411	CPR Prompt® Rescue and
	NG Tube & Trach Skills	L1002040	Practice Aid
LFU11/4U	ING TUDE & TRUCH SKIIIS		Tractice Ala

Nasco Fort Atkinson

901 Janesville Avenue, P.O. Box 901 Fort Atkinson, Wisconsin 53538-0901 1-800-558-9595

eNasco.com • E-mail: lifeform@eNasco.com COPYRIGHT © NASCO 2001 PRINTED IN U.S.A. NP 62-83/RV 3-12